



LEAP PAL PARTS + CONSUMABLES

PEEK TUBING DISCOUNT PROMO - 40% OFF



PEEK TUBING PRODUCT LISTING

Product Code	Description	List Price	Discount	Promo Price
VAL.JR-T-5804-100	Solid Color Coded PEEK Tubing (Black) - 1/16" OD x .004" ID x 100FT	\$435.50	\$174.20	\$261.30
VAL.JR-T-5804ST-100	Striped Color Coded PEEK Tubing (Black) - 1/16" OD x .004" ID x 100FT	\$435.50	\$174.20	\$261.30
VAL.JR-T-5999-100	Striped Color Coded PEEK Tubing (Red) - 1/16" OD x .005" ID x 100FT; 6700 PSI	\$390.00	\$156.00	\$234.00
VAL.JR-T-6000-100	Striped Color Coded PEEK Tubing (Yellow) - 1/16" OD x .007" ID x 100FT; 6400 PSI	\$325.00	\$130.00	\$195.00
VAL.JR-T-6001-100	Striped Color Coded PEEK Tubing (Blue) - 1/16" OD x .010" ID x 100FT; 6100 PSI	\$325.00	\$130.00	\$195.00
VAL.JR-T-6002-100	Striped Color Coded PEEK Tubing (Orange) - 1/16" OD x .020" ID x 100FT; 5000 PSI	\$325.00	\$130.00	\$195.00
VAL.JR-T-6003-100	Striped Color Coded PEEK Tubing (Green) - 1/16" OD x .030" ID x 100FT; 3800 PSI	\$325.00	\$130.00	\$195.00
VAL.JR-T-60031-100	Striped Color Coded PEEK Tubing (Grey) - 1/16" OD x .040" ID x 100FT; 2600 PSI	\$325.00	\$130.00	\$195.00
VAL.JR-T-6007-100	Solid Layer Color Coded PEEK Tubing (Red) - 1/16" OD x .005" ID x 100FT; 6700 PSI, 460 BAR	\$390.00	\$156.00	\$234.00
VAL.JR-T-6008-100	Solid Layer Color Coded PEEK Tubing (Yellow) - 1/16" OD x .007" ID x 100FT; 6400 PSI, 440 BAR	\$325.00	\$130.00	\$195.00
VAL.JR-T-6009-100	Solid Layer Color Coded PEEK Tubing (Blue) - 1/16" OD x .010" ID x 100FT; 6100 PSI, 420 BAR	\$325.00	\$130.00	\$195.00
VAL.JR-T-6010-100	Solid Layer Color Coded PEEK Tubing (Orange) - 1/16" OD x .020" ID x 100FT; 5000 PSI, 345 BAR	\$325.00	\$130.00	\$195.00
VAL.JR-T-6010-500	Tubing, PEEK, 1/16"OD x 0.020", Orange 500 ft Spool	\$1,444.00	\$577.60	\$866.40
VAL.JR-T-6011-100	Solid Layer Color Coded PEEK Tubing (Green) - 1/16" OD x .030" ID x 100FT; 3800PSI, 262 BAR	\$325.00	\$130.00	\$195.00

LEAP PAL PARTS + CONSUMABLES

Product Code	Description	List Price	Discount	Promo Price
VAL.JR-T-62015-100	Solid Layer Color Coded PEEK Tubing (Grey) - 1/16" OD x .015" ID x 100FT	\$325.00	\$130.00	\$195.00
VAL.JR-TD-5804-100	Dual Layer Color Coded PEEK Tubing (Black) - 1/16" OD x .004" ID x 100FT; 6700 PSI	\$435.50	\$174.20	\$261.30
VAL.JR-TD-6007-100	Dual Layer Color Coded PEEK Tubing (Red) - 1/16" OD x .005" ID x 100FT; 6700 PSI	\$317.20	\$126.88	\$190.32
VAL.JR-TD-6008-100	Dual Layer Color Coded PEEK Tubing (Yellow) - 1/16" OD x .007" ID x 100FT; 6400 PSI	\$317.20	\$126.88	\$190.32
VAL.JR-TD-6009-100	Dual Layer Color Coded PEEK Tubing (Blue) - 1/16" OD x .010" ID x 100FT; 6100 PSI	\$317.20	\$126.88	\$190.32
VAL.JR-TD-6010-100	Dual Layer Color Coded PEEK Tubing (Orange) - 1/16" OD x .020" ID x 100FT; 5000 PSI	\$317.20	\$126.88	\$190.32
VAL.JR-TD-6011-100	Dual Layer Color Coded PEEK Tubing (Green) - 1/16" OD x .030" ID x 100FT; 3800 PSI	\$317.20	\$126.88	\$190.32
VAL.TPK.102-100	Natural PEEK Tubing - 360 µm OD x .002" ID x 100FT; 5200 PSI	\$585.00	\$234.00	\$351.00
VAL.TPK.103-100	Natural PEEK Tubing - 360 µm OD x .003" ID x 100FT; 5000 PSI	\$585.00	\$234.00	\$351.00
VAL.TPK.104-100	Natural PEEK Tubing - 360 µm OD x .004" ID x 100FT; 4400 PSI	\$585.00	\$234.00	\$351.00
VAL.TPK.105-100	Natural PEEK Tubing - 360 µm OD x .005" ID x 100FT; 4000 PSI	\$455.00	\$182.00	\$273.00
VAL.TPK.106-100	Natural PEEK Tubing - 360 µm OD x .006" ID x 100FT; 3500 PSI	\$585.00	\$234.00	\$351.00
VAL.TPK.502-100	Natural PEEK Tubing - 1/32" OD x .0025" ID x 100FT; 6600 PSI	\$390.00	\$156.00	\$234.00
VAL.TPK.505-100	Natural PEEK Tubing - 1/32" OD x .005" ID x 100FT; 6000 PSI	\$390.00	\$156.00	\$234.00
VAL.TPK.507-100	Natural PEEK Tubing - 1/32" OD x .007" ID x 100FT; 5600 PSI	\$403.00	\$161.20	\$241.80
VAL.TPK.510-100	Natural PEEK Tubing - 1/32" OD x .010" ID x 100FT; 4900 PSI	\$338.00	\$135.20	\$202.80
VAL.TPK.515-100	Natural PEEK Tubing - 1/32" OD x .015" ID x 100FT; 3900 PSI	\$338.00	\$135.20	\$202.80
VAL.TPK102-100	Natural PEEK Tubing - 1/16" OD x .002" ID x 100FT; 7000 PSI	\$455.00	\$182.00	\$273.00
VAL.TPK103-100	Natural PEEK Tubing - 1/16" OD x .003" ID x 100FT; 6800 PSI	\$455.00	\$182.00	\$273.00
VAL.TPK104-100	Natural PEEK Tubing - 1/16" OD x .004" ID x 100FT; 6700 PSI	\$455.00	\$182.00	\$273.00
VAL.TPK105-100	Natural PEEK Tubing - 1/16" OD x .005" ID x 100FT; 6700 PSI	\$455.00	\$182.00	\$273.00
VAL.TPK106-100	Natural PEEK Tubing - 1/16" OD x .006" ID x 100FT; 6500 PSI	\$390.00	\$156.00	\$234.00
VAL.TPK107-100	Natural PEEK Tubing - 1/16" OD x .007" ID x 100FT; 6400 PSI	\$390.00	\$156.00	\$234.00
VAL.TPK110-100	Natural PEEK Tubing - 1/16" OD x .010" ID x 100FT; 6100 PSI	\$312.00	\$124.80	\$187.20
VAL.TPK115-100	Natural PEEK Tubing - 1/16" OD x .015" ID x 100FT; 5500 PSI	\$364.00	\$145.60	\$218.40
VAL.TPK120-100	Natural PEEK Tubing - 1/16" OD x .020" ID x 100FT; 5000 PSI	\$312.00	\$124.80	\$187.20
VAL.TPK130-100	Natural PEEK Tubing - 1/16" OD x .030" ID x 100FT; 3800 PSI	\$312.00	\$124.80	\$187.20
VAL.TPK140-100	Natural PEEK Tubing - 1/16" OD x .040" ID x 100FT; 2600 PSI	\$364.00	\$145.60	\$218.40
VAL.TPK155-100	Natural PEEK Tubing - 1/16" OD x .055" ID x 100FT; 800 PSI	\$312.00	\$124.80	\$187.20
VAL.TPK260-100	Natural PEEK Tubing - 1/8" OD x .060" ID x 100FT; 3600 PSI	\$520.00	\$208.00	\$312.00