



### General Arrangement Drawings Deepwell Plates for Chromatography-Round Well

ITEMCODE	Description	A	B	C	D	E	F	G1	G2	H	J	K	L	M	N
LPP.1300RRPP	96-Well MicroWell Microplate, round well-2mL	127.8	85.5	125.4	84.3	2.5	31.6	Ø8.4	Ø7.2	14.4	11.2	9	29.1	1.4	2.5
LPP.2000RUPP	96-Well MicroWell Microplate, round well-1.3mL	127.8	85.5	123.6	81.6	2.5	43.6	Ø8.5	Ø7.2	14.3	11.3	9	41.5	1.4	2.5

#### LPP.1300RRPP/LPP.2000RUPP



**Notes.**

Drawings are for instrument set-up purposes and are not manufacturing drawings

LEAP PAL Parts is providing information that has been supplied by the original manufacturer and is not responsible for subsequent changes or modifications.