


General Arrangement Drawings Microplates for Chromatography-Certified

ITEMcode	Description	A	B	F	G1	H	J	K	L	M	N
		BASE LENGTH	BASE WIDTH	OVERALL HEIGHT	DIA TOP	A - I LOCATION (X)	A - I LOCATION (Y)	WELL TO WELL SPACING	WELL DEPTH	BOTTOM OF WELL DISTANCE	FLANGE SHORT SIDE
LPP.350RFPF	Round Well-Flat bottom	128.00	85.40	14.70	7.05	14.70	11.50	9.00	13.90	3.00	3.00
LPP.270RUPP	Round Well-Round bottom	128.00	85.40	14.70	7.05	14.70	11.50	9.00	13.90	3.00	3.00
LPP.220RVPP	Round Well-V bottom	128.00	85.40	14.70	7.05	14.70	11.50	9.00	13.90	3.00	3.00



Notes.

Drawings are for instrument set-up purposes and are not manufacturing drawings

LEAP PAL Parts is providing information that has been supplied by the original manufacturer and is not responsible for subsequent changes or modifications.

All dimensions mm