

NUNC NALGENE PLATES

AVAILABLE NOW AT COMPETITIVE PRICES

The highest quality, most durable microplates in the industry providing superior performance, long product life and dependably low extractables.

FEATURES

- Manufactured in Class 8 clean room, free of contamination
- Made with virgin, medical grade, polypropylene resin
- All plates meet ANSI/SLAS standards

BENEFITS

- Highest sterility assurance in the industry
- Low binding characteristics, allow complete sample recovery
- Optimized for robotics, readers and liquid handling systems

Plates

ITEM CODE

DESCRIPTION

•LPP.450RVPP	-NUNC MicroWell Plate with 96 Round wells, 8mm, Conical Bottom Shape, Non-sterile, 400 uL Working Volume, Raised Rim, 120-Pack
•LPP.500RUPP	-NUNC MicroWell Plate with 96 Round wells, 8 mm, Round Bottom Shape, Non-sterile, 450 uL Working Volume, Raised Rim, 120-Pack
•LPP.1300RRPP	-NUNC MicroWell Plate with 96 Round wells, 8 mm, Round Bottom Shape, Non-sterile, 1000 uL Working Volume, Raised Rim, 50-Pack
•LPP.2000RUPP	-NUNC DeepWell Plate with 96 Round wells, 8 mm, Round Bottom Shape, Non-sterile, 1900 uL Working Volume, Raised Rim, 60-Pack
•LPP.2000SVPP	-NUNC Plate with 96 Square wells, Conical Bottom Shape, Non-sterile, 1900 uL Working Volume, Raised Rim, 50-Pack

***For information on available discounts, please speak with a member of the sales team.**

Well Shape




Round

- Wide range of sealing options
- Ideal for sample mixing
- Ideal for samples requiring vortexing



Square

- Maximizes volume
- Easier for automated plate changes



Well Bottom



Round

- Ideal for mixing




Conical

- Improves sample recovery and decreased dead volume



Flat

- Useful for microplate readers and absorbance readings



Cap Mats

ITEM CODE

DESCRIPTION

•LPP.209112V	-L-MARK CapMat 96 well, Round, Silicone/PTFE coated, 7mm, pre-slit, 5-Pack
•LPP.212512V	-L-MARK CapMat 96 well, Round, Silicone/PTFE coated, 8mm, pre-slit, 5-Pack
•LPP.96R8MMCLSPS	-NUNC CapMat 96 well, Round, Silicone, 8mm, pre-slit, 10-Pack
•POR.500179	-Porvair CapMat, 96 Well, Round, Silicone, 8mm, cross-score, 50-Pack
•LPP.208922V	-L-MARK CapMat 96 well, Square, Silicone/PTFE coated, pre-slit, 5-Pack

CALL FOR SAMPLES TODAY!

LEAP PAL PARTS + CONSUMABLES®
(877) 725 1007 | palparts.com

LEAP
PAL
PARTS
+ CONSUMABLES