

# Connections | GC Ferrules SilTite® FingerTite

The Smart Alternative. One Nut, One Ferrule, All Systems.

SilTite® FingerTite is a GC ferrule system that has been designed to be installed only with the force of you fingers, delivering an easy, leak and maintenance free seal.

### Features and benefits

- Simplified installation process saves time.
- Leak free no column creep.
- Cannot overtighten no broken columns.
- Chemically inert no contamination.
- Metal ferrule no temperature limit.

# Recommended applications

MS, high throughput laboratories.

## Product specifications

Kits available for Agilent, PerkinElmer, Shimadzu, Thermo Scientific, and Varian/Bruker GC-MS systems.

Ferrules available with 0.4, 0.5, and 0.7 mm hole sizes. Blanking ferrules, replacement female nuts, and replacement Agilent base seals also available.

For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com